


putting PPE **ON**



(ACE)



Contact gown



N95 Respirator



Face shield or goggles



Gloves

## What You Need

### Gown + Gloves

1 **Remove** any personal items and jewelry and put in secure location, not in pockets.

2 **Sanitize** hands.

3 **Put on** contact gown outside room.

**Gown** opens in back.

**Tie** the back of the gown.

4 **Put on** gloves over the cuffs of the gown.

### Mask

5 **Put on** N95, ensuring proper seal.

**Ensure** straps are not crossing.

6 **Place** hands over the front of the N95. **Breathe** an easy deep breath in and out.

If you **feel** air escape the edges, **refit** and **repeat**.

NO NO NO NO NO NO

YES

7 **Sanitize** gloves.

### Face Shield

8 **Put on** face shield.

9 **Sanitize** gloves.

10 **ENTER** room

**DO NOT** enter the room if you do not achieve a proper respirator seal.

11 **Do not** touch face or re-adjust N95 or face shield inside room.